


Bacterial Transformation  
Mohana Greeshma Gudhimalla  
04/14/2023.

The experimental results were successful in all plates except in the plate where there is no pGLO, LB, and Amp. The results predicted were 0 because the ampicillin antibiotic would not provide an environment for growth without the presence of pGLO, but the experiment resulted in colonies too numerous to count (TNTC) because of non-functional ampicillin.

Plates Characteristics	Plasmid added	Colony number	Colony fluorescence
- , LB ,	No	TNTC	no
- , LB , Amp	No	TNTC	no
+ , LB , Amp	Yes	TNTC	no
+ , LB , Amp , Ara	Yes	2	Yes



- DNA plasmid concentration: 0.05  $\mu\text{g}/\mu\text{l}$
- Cells incubated with 400  $\mu\text{l}$   $\text{CaCl}_2$  transformation solution and split into two tubes
- 5  $\mu\text{l}$  pGLO plasmid solution to the tube of +pGLO
- 200  $\mu\text{l}$  LB broth
- 100  $\mu\text{l}$  cells spread on agar
- 2 colonies of transformants on the LB/amp plate

Total amount of pGLO DNA used (0.5 pts) = 0.25  $\mu\text{g}$

Fraction of DNA used (0.5 pts) = 0.2469

Micrograms of DNA spread on the plates (0.5 pts) = 0.06173  $\mu\text{g}$

Transformation efficiency (0.5 pts) = 32.4 transformants/ $\mu\text{g}$